

Anti-MTPN antibody (1-118 aa) (STJ11105943)

STJ11105943

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	136319
Gene Symbol	MTPN
Uniprot ID	MTPN_HUMAN
Immunogen	
Immunogen Region	1-118 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human MTPN (NP_665807.1).
Immunogen Sequence	MCDKEFMWALKNGDLDEVKD YVAKGEDVNRTEGGRKPLH YAADCGQLEILEFLLKGAD INAPDKHHITPLLSAVYEGH VSCVKLLLSKGADKTVKGPD GLTAFEATDNQAIKALLQ