

Anti-RELL2 antibody (44-303 aa) (STJ11105950)

STJ11105950

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	285613
Gene Symbol	RELL2
Uniprot ID	RELL2_HUMAN
Immunogen	
Immunogen Region	44-303 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 44-303 of human RELL2 (NP_776189.3).
Immunogen Sequence	YRCRTSRGSEPDDAQLQPPE DDDMNEDTVERIVRCIIQNE ANAEALKEMLGDSEGETVQ LSSVDATSSLQDGAPSHHHT VHLGSAAPCLHCSRKRPL VRQGRSKEGKSRPRTGETTV FSVGRFRVTHIEKRYGLHEH RDGSPTDRSWGSGGQDPGG GQGSGGGQPKAGMPAMERLP PERPQPQLASPPVQNGGLR DSSLTPRALEGNPRASAEPT LRAGGRGPSPGLPTQEANG