

Anti-WDR38 antibody (1-314 aa) (STJ11105951)

STJ11105951

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	401551
Gene Symbol	WDR38
Uniprot ID	WDR38_HUMAN
Immunogen	
Immunogen Region	1-314 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-314 of human WDR38 (NP_001038941.1).
Immunogen Sequence	MNSGVPATLAVRRVKFFGQH GGEVNSSAFSPDGQMLLTGS EDGCVYGWETRSGQLLWRLG GHTGPVKFCRFSPDGHLFAS ASCDCTVRLWDVARAKCLRV LKGHQRSVETVSFSPDSRQL ASGGWDKRVMLWDVQSGQML RLLVGHDRDSIQSSDFSPSTVN CLATGSWDSTVHIWDLRMVT PAVSHQALEGHSANISCLCY SASGLLASGSDKTIHIWKP TTSSLLIQLKGHVTWVKSI