

Anti-PAPSS2 antibody (306-487 aa) (STJ11105952)

STJ11105952

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9060
Gene Symbol	PAPSS2
Uniprot ID	PAPS2_HUMAN
Immunogen	
Immunogen Region	306-487 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 306-487 of human PAPSS2 (NP_004661.2).
Immunogen Sequence	DKTRLEGCSKFVLAHGGRRV AILRDAEFYEHKEERCSR V WGTCTKHPIKMMVMSGDW LVGGDLQVLEKIRWNDGLDQ YRLTPLELKQCKEMNADAV FAFQLRNPVHNGHALLMQDT RRRLLERGYKHPVLLHPLG GWTKDDDVPLDWRMKQHAHV LEEGVLPKSTIVAIFPSPM LY