

Anti-KERA antibody (21-352 aa) (STJ11105954)

STJ11105954

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	11081
Gene Symbol	KERA
Uniprot ID	KERA_HUMAN
Immunogen	
Immunogen Region	21-352 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 21-352 of human KERA (NP_008966.1).
Immunogen Sequence	RSVRQVYEVHDSDDWTIHDF ECPMECFCPPSFPTALYCEN RGLKEIPAIPSRWYLYLQN NLIETIPEKPFENATQLRWI NLNKNKITNYGIEKGALSQL KLLFLFLEDNELEEVPSPL PRSLEQLQLARNKVSRIQQG TFSNLENLTLLDLQNNKLVD NAFQRDTFKGLKNLMQLNMA KNALRNMPPRLPANTMQLFL DNNSIEGIPENYFNVIPKVA FLRLNHNKLSDEGLPSRGF