

Anti-B9D2 antibody (1-175 aa) (STJ11105955)

STJ11105955

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	80776
Gene Symbol	B9D2
Uniprot ID	B9D2_HUMAN
Immunogen	
Immunogen Region	1-175 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human B9D2 (NP_085055.2).
Immunogen Sequence	MAEVHVGIIIGASGFESS LFCKWGIHTGAAWKLLSGVR EGQTQVDTPQIGDMAYWSHP IDLHFATKGLQGWPRLHFQV WSQDSFGRCQLAGYGFCVHP SSPGTHQLACPTWRPLGSWR EQLARAFVGGGPQLLHGDTI YSGADRYRLHTAAGGTVHLE IGLLLRNFDTRYGVEC