

Anti-TTLL7 antibody (1-260 aa) (STJ11105957)

STJ11105957

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	79739
Gene Symbol	TTLL7
Uniprot ID	TTLL7_HUMAN
Immunogen	
Immunogen Region	1-260 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TTLL7 (NP_078962.4).
Immunogen Sequence	MPSLPQEGVIQGPSPLDLNT ELPYQSTMKRKVRKKKKKGT ITANVAGTKFEIVRLVIDEM GFMKTPDEDETSNLIWCDSA VQQEKISELQNYQRINHFPG MGEICRKDFLARNTMKMIKS RPLDYTFVPRTWIFPAEYTD FQNYVKELKKRKQKTFIVK PANGAMGHGISLIRNGDKLP SQDHLIVQEYIEKPFLMEGY KFDLRIYILVTSCDPLKIFL YHDGLVRMGTEKYIPPNES