

Anti-GNG11 antibody (1-73 aa) (STJ11105959)

STJ11105959

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	2791
Gene Symbol	GNG11
Uniprot ID	GBG11_HUMAN
Immunogen	
Immunogen Region	1-73 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-73 of human GNG11 (NP_004117.1).
Immunogen Sequence	MPALHIEDLPEKEKLKMEVE QLRKEVKLQRQQVSKCSEEI KNYIEERSGEDPLVKGIPED KNPFKKEKGSCVIS