

Anti-TRIM10 antibody (136-395 aa) (STJ11105964)
STJ11105964

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10107
Gene Symbol	TRIM10
Uniprot ID	TRI10_HUMAN
Immunogen	
Immunogen Region	136-395 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 136-395 of human TRIM10 (NP_439893.2).
Immunogen Sequence	EDAAAPYREQIHKCLKLRK EREEIQEIQSRENKRMQVLL TQVSTKRQQVISEFAHLRKF LEEQQSILLAQLESQDGDIL RQRDEFDLLVAGEICRFSAL IEELEEKNERPARELLTDIR STLIRCETRKCRKPVAVSPE LGQRIRDFPQQALPLQREMK MFLEKLCFELDYEPAHISLD PQTSHPKLLLSSEDHQRAQFS YKWQNSPDNPQRFRATCVL AHTGITGGRHTWWMARVP