

Anti-TSC22D1 antibody (1-144 aa) (STJ11105967)

STJ11105967

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	8848
Gene Symbol	TSC22D1
Uniprot ID	T22D1_HUMAN
Immunogen	
Immunogen Region	1-144 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human TSC22D1 (NP_006013.1).
Immunogen Sequence	MKSQWCRPVAMD LGVYQLRH FSISFLSLLGTENASVRLD NSSSGASVVAIDNKIEQAMD LVKSHLMYAVREEVEVLKEQ IKELIEKNSQLEQENLLKT LASPEQLAQFQAQLQTGSPPTTQPQGTTPPAQPASQGS GPTA