

Anti-PTPN9 antibody (294-593 aa) (STJ11105969)
STJ11105969

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5780
Gene Symbol	PTPN9
Uniprot ID	PTN9_HUMAN
Immunogen	
Immunogen Region	294-593 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 294-593 of human PTPN9 (NP_002824.1).
Immunogen Sequence	YVNARQKQGIYEEYEDIRRE NPVGTFHCMSPGNLEKNRY GDVPCLDQTRVKLTKRSGHT QTDYINASFMDGYKQKNAYI GTQGPLENTYRDFWLMVWEQ KVLVIVMTTRFEEGGRRKCG QYWPLEKDSRIIRFGFLTVTN LGVENMNHYYKKTLEIHNTE ERQKRQVTHFQFLSWPDYGV PSSAASLIDFLRVVRNQQL AVSNMGARSKGQCPEPPIV HCSAGIGRTGTFCSLDIDL