

Anti-ARHGEF28 antibody (1050-1250 aa) (STJ11105970)

STJ11105970

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64283
Gene Symbol	ARHGEF28
Uniprot ID	ARG28_HUMAN
Immunogen	
Immunogen Region	1050-1250 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1050-1250 of human ARHGEF28 (NP_001073948.2).
Immunogen Sequence	KNQKWLEILNKIENKTYTKL KNGHVFRKQALMSEERTLLY DGLVYWKATGRFKDILALL LTDVLLFLQEKDQKYIFAAV DQKPSVISLQKLIAREVANE ERGMFLISASSAGPEMYEIH TNSKEERNNWMRRIQQAVES CPEEKGGRTSESDKRAE ARVAKIQCCQEILTNDQDQI CAYLEELHIYAEELGELSGF E



Western blot analysis of lysates from HeLa cells, using ARHGEF28 Rabbit pAb (STJ11105970) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
 Lysates/proteins: 25 μ g per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081