

Anti-ADCY4 antibody (808-1077 aa) (STJ11105971)

STJ11105971

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	196883
Gene Symbol	ADCY4
Uniprot ID	ADCY4_HUMAN
Immunogen	
Immunogen Region	808-1077 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 808-1077 of human ADCY4 (NP_640340.2).
Immunogen Sequence	RQNEYCYRLDFLWKKLRQE REETETMENLTRLLENVLP AHVAPQFIGQNRNEDLYHQ SYECVCLFASVPDFKEFYS ESNINHEGLECLRLLENEIA DFDELLSKPKFSGVEKIKTI GSTYMAATGLNATSGQDAQQ DAERSCSHLGTMVEFAVALG SKLDVINKHSFNNFRLRVGL NHGPVWAGVIGAQPQYDIW GNTVNVASRMESTGVLGKIQ VTEETAWALQSLGYTCYSR