

Anti-CALML3 antibody (1-149 aa) (STJ11105972)

STJ11105972

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	810
Gene Symbol	CALML3
Uniprot ID	CALL3_HUMAN
Immunogen	
Immunogen Region	1-149 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human CALML3 (NP_005176.1).
Immunogen Sequence	MADQLTEEQVTEFKEAFSLF DKDGDGCITTRELGTVMRSL GQNPTEAELRDMMSSEIDRDG NGTVDFPEFLGMMARKMKDT DNEEEIREAFRVFDKDGNGF VSAAELRHVMTRLGEKLSDE EVDEMIRAADTDGGQVNYE EFVRVLVSK