

Anti-GCM2 antibody (287-506 aa) (STJ11105978)

STJ11105978

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9247
Gene Symbol	GCM2
Uniprot ID	GCM2_HUMAN
Immunogen	
Immunogen Region	287-506 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 287-506 of human GCM2 (NP_004743.1).
Immunogen Sequence	PYPTLYKDYKSTIPNDTDWVH LNTLQCNVNSYSSYERSFDF TNKQHGWKPALGKPSLVERT NHGQFQAMATRPYYNPELPC RYLTTPPPGAPALQTVITTT TKVSYQAYQPPAMKYSDSVR EVKSLSSCNYPEDTGMSVY PEPWGPPVTVTRAASPSGPP PMKIAGDCRAIRPTVAIPHE PVSSRTDEAETWDVCLSLG SAVSYSDRVGPFFTYNNEDF