

Anti-GSPT2 antibody (1-190 aa) (STJ11105982)

STJ11105982

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23708
Gene Symbol	GSPT2
Uniprot ID	ERF3B_HUMAN
Immunogen	
Immunogen Region	1-190 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human GSPT2 (NP_060564.2).
Immunogen Sequence	MDSGSSSSDSAPDCWDQVDM ESPGSAPSGDGVSSAVAEQ REPLSSAFSRKLVNNAKPFV PNVHAAEFVPSFLRGPTQPP TLPAGSGSNDCTGAGYPQ GKRMGRGAPVEPSREEPLVS LEGSNSAVTMELSEPWENG EVEMALEESWEHSKEVSEAE PGGSSGDSGPPEESGQEMM EEEKEIRKSK