

Anti-KIF1A antibody (1596-1690 aa) (STJ11105984)

STJ11105984

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	547
Gene Symbol	KIF1A
Uniprot ID	KIF1A_HUMAN
Immunogen	
Immunogen Region	1596-1690 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1596-1690 of human KIF1A (NP_004312.2).
Immunogen Sequence	FVVRRPYAYMYNSDKDTVE RFLNLATAQVEYSEDQQAM LKTPNTFAVCTEHRGILLQA ASDKMDHDLWLYAFNPPLAGT IRSKLSRRRSQMRV

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081