

Anti-SLC24A1 antibody (330-452 aa) (STJ11105987)

STJ11105987

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9187
Gene Symbol	SLC24A1
Uniprot ID	NCKX1_HUMAN
Immunogen	
Immunogen Region	330-452 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 330-452 of human SLC24A1 (O60721).
Immunogen Sequence	GSSPAETKAFTAAWSLRNPS PRTSVSAIKTAPAIWRLAK KPSTAPSTSTPTVRAKLTM QVHHCVVVKPTPAMLTTPSP SLTTALLPEELSPSPSVLPP SLPDLHPKGEYPPDLFSVEE RRQ