

Anti-HTR6 antibody (321-440 aa) (STJ11105988)

STJ11105988

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3362
Gene Symbol	HTR6
Uniprot ID	5HT6R_HUMAN
Immunogen	
Immunogen Region	321-440 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 321-440 of human HTR6 (NP_000862.1).
Immunogen Sequence	PLFMRFKRALGRFLPCPRC PRERQASLASPSLRTSHSGP RPGLSLQQVLPLPLPPSDS DSDAGSGGSSGLRLTAQLLL PGEATQDPPLPTRAAAVNF FNIDPAEPELRPHPLGIPTN