

Anti-CCL23 antibody (22-120 aa) (STJ11105997)

STJ11105997

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6368
Gene Symbol	CCL23
Uniprot ID	CCL23_HUMAN
Immunogen	
Immunogen Region	22-120 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-120 of human CCL23 (NP_665905.1).
Immunogen Sequence	RVTKDAETEFMMSKLPENP VLLDRFHATSADCCISYTPR SIPCSLLESYFETNSECSKP GVIFLTKKGRRFCANPSDKQ VQVCVRMLKLDTRIKTRKN

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081