

Anti-JPH3 antibody (500-727 aa) (STJ11106001)

STJ11106001

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	57338
Gene Symbol	JPH3
Uniprot ID	JPH3_HUMAN
Immunogen	
Immunogen Region	500-727 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 500-727 of human JPH3 (NP_065706.2).
Immunogen Sequence	HFSRQVSVDEERGDDIQLL EGRAGDCARSSWGEEQAGGS RGVRS GALRGLLVDDFRTR GSGRKQPGNPKPRERRTESP PVFTWTS HHASNHSPGGSRLLELQEEKLSNYRMEMKPLL RMETHPQKRRYSKGGACRGL GDDHRPEDRGFGVQRLRSKA QNKENFRPASSAEPVQKLA SLRLGGAEPRLLRWDLTFS P QKSLPVALESDEENGDELK SSTGSAPI