

Anti-HEPH antibody (300-580 aa) (STJ11106008)

STJ11106008

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9843
Gene Symbol	HEPH
Uniprot ID	HEPH_HUMAN
Immunogen	
Immunogen Region	300-580 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 300-580 of human HEPH (NP_620074.3).
Immunogen Sequence	ATYCSDPASVDKETFOES NRMHAINGFVFGNLPENMC AQKRVAVHFLFGMGNEIDVHT AFFHGQMLTTRGHHTDVANI FPATFVTAEMVPWEPGTWLI SCQVNSHFRDGMQALYKVKS CSMAPPVDLLTGKVRQYFIE AHEIQWDYGPMSGHDGSGKN LREPGSISDKFFQKSSSRIG GTYWKVRYEAFQDETQEKM HLEEDRHLGILGPVIRAEVG DTIQVVFYNRASQPFSMQP