

Anti-KCNJ2 antibody (50-150 aa) (STJ11106018)

STJ11106018

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3759
Gene Symbol	KCNJ2
Uniprot ID	KCNJ2_HUMAN
Immunogen	
Immunogen Region	50-150 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human KCNJ2 (NP_000882.1).
Immunogen Sequence	KDGHGCVQFINVGEKGQRYL ADIFTTCVDIRWRWMLVIFC LAFVLSWLFFGCVFWLIALL HGDLDASKEGKACVSEVNSF TAAFLFSIETQTTIGYGFRC V