

Anti-CLDN12 antibody (1-100 aa) (STJ11106019)

STJ11106019

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9069
Gene Symbol	CLDN12
Uniprot ID	CLD12_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CLDN12 (NP_001172001.1).
Immunogen Sequence	MGCRDVHAATVLSFLCGIAS VAGLFAGTLLPNWRKLRIT FNRNEKNLTVYTGLWVKCAR YDGSSDCLMYDTTWYSSVDQ LDLRVLQFALPLSMLIAMGA