

Anti-MT-ND4L antibody (1-98 aa) (STJ11106020)

STJ11106020

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	17720
Gene Symbol	Mtnd4l
Uniprot ID	NU4LM_MOUSE
Immunogen	
Immunogen Region	1-98 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-98 of mouse MT-ND4L (NP_904336.1).
Immunogen Sequence	MPSTFFNLTMAFSLSLGTL MFRSHLMSTLLCLEGMVLSL FIMTSVTSLSNSMSMSSMPIP ITILVFAACEAAVGLALLVK VSNTYGTDYVQNLNLLQC

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081