

Anti-UBN1 antibody (1000-1100 aa) (STJ11106021)

STJ11106021

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	170644
Gene Symbol	Ubn1
Uniprot ID	UBN1_MOUSE
Immunogen	
Immunogen Region	1000-1100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1000-1100 of mouse UBN1 (NP_080942.1).
Immunogen Sequence	SVTSSTSLPKGTGGAVLLSN TSSL SLLSSSYKSNNPKLPG AMNSNSLGIITQFPLHVLSF NADSSAKAGVSKDAIVTGPA PGTFHHGLSHSLLAGLHSSP P