

Anti-C15orf43 antibody (1-218 aa) (STJ11106031)

STJ11106031

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	74401
Gene Symbol	Terb2
Uniprot ID	TERB2_MOUSE
Immunogen	
Immunogen Region	1-218 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of mouse C15orf43 (NP_083190.1).
Immunogen Sequence	MFQGGQRGWFCGSVSQDLRQI WEDEGGMVSDVKAADFLFSC DASHPDLRIYQSLEYIEDN ATVFHAYYLAAIANTEMKNS VALGHFVLPPACQLQKEIRRK IGSFIWEQDEKFQIEKHDRM ASSDKENIRPTPEHKQELSK SAEHHLTRTPVIEKQMCPL HSYPVNNMVTGYISIDALEK FLGELHDFTPGSSGYLAYHI QDEINMSAIKNKLRRKLS