

Anti-FGL1 antibody (STJ11106037)

STJ11106037

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	2267
Gene Symbol	FGL1
Uniprot ID	FGL1_HUMAN
Immunogen	
Immunogen Region	
Specificity	A synthetic Peptide of human FGL1
Immunogen Sequence	



Western blot analysis of lysates from fetal liver, using FGL1 Rabbit pAb (STJ11106037).
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000656) at 1:10000 dilution.
Lysates/proteins: 25 µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081