

Anti-Beta III-Tubulin antibody (40-260 aa) (STJ11106048)
STJ11106048

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/FC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:100-1:500
Range	IF/CC:1:50-1:200 FC:1:20-1:50 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10381
Gene Symbol	TUBB3
Uniprot ID	TBB3_HUMAN
Immunogen	
Immunogen Region	40-260 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 40-260 of human Beta III-Tubulin/Beta 3-Tubulin (NP_006077.2).
Immunogen Sequence	SDLQLERISVYYNEASSHIKY VPRAILVDLEPGTMDSVRSQ AFGHLFRPDNFIQSGAGN NWAKGHYTEGAELVDSVLDV VRKECENCDCQLQGFLTHSL GGGTGSGMGTLISKVREEY PDRIMNTFSVVPSPKVSDTV VEPYNATLSIHQLVENTDET YCIDNEALYDICFRTLKLAT PTYGDLNHLVSATMSGVTTS LRFPGQLNADLRKLAVNMVP F



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081