

Anti-URI1 antibody (231-531 aa) (STJ11106058)

STJ11106058

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	19777
Gene Symbol	Uri1
Uniprot ID	RMP_MOUSE
Immunogen	
Immunogen Region	231-531 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 231-531 of mouse URI1 (NP_035404.4).
Immunogen Sequence	TVIANGEDRVSEEEKEGAD TGVNVVSPVTDSSAASSCKR RAGNAGLPNGQVNSLNYSVN GSNSYHSNKDDDEEEEDDDD DDDEDDNESDHAISADNSI PTIYFSHTVEPKRVRINTGK NTTLKFSEKKEAKRKRKSG AGSHATHELPAIKSPADIYR VFVDVVNGEYVPRKSILKSR SRENSVCSDTSESSAADVED RRLLRSTSSSEEAVATEAGG SSLDELQENHPKKPLPSGV



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081