

Anti-PPP1R15B antibody (494-713 aa) (STJ11106060)

STJ11106060

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84919
Gene Symbol	PPP1R15B
Uniprot ID	PR15B_HUMAN
Immunogen	
Immunogen Region	494-713 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 494-713 of human PPP1R15B (NP_116222.4).
Immunogen Sequence	ATIQTAAARIVPEEPSDSEKD LSGKSDLENSSQSGSLPETP EHSSEEDDWESSADEAESL KLWNSFCNSDDPYNPLNFKA PFQTSGENEGCRDSKTPSE SIVAISECHTLLSCKVQLLG SQESECPSVQRDVLSSGGRH THVKRKKVTFLEEVTEYYIS GDEDRKGPWEEFARDGCRFQ KRIQETEDAIGYCLTFEHRE RMFNRLQGTCTFKGLNVLKQC