

Anti-FOXO1 antibody (184-373 aa) (STJ11106069)

STJ11106069

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2304
Gene Symbol	FOXO1
Uniprot ID	FOXO1_HUMAN
Immunogen	
Immunogen Region	184-373 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 184-373 of human FOXO1 (NP_004464.2).
Immunogen Sequence	AVPAARPPYPGAVYAGYAPP SLAAPPVYYPAAASPGPCR V FGLVPERPLSPELGPAPSGP GGSCAFASAGAPATTGYQP AGCTGARPNPSAYAAAYAG PDGAYPQGAGSAIFAAAGRL AGPASPPAGGSSGGVETTVD FYGRTSPGQFGALGACYNPG QLGGASAGAYHARHAAAYP GGIDRFVSAM