

Anti-CDY2B antibody (1-280 aa) (STJ11106070)

STJ11106070

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	203611/9426
Gene Symbol	CDY2A.CDY2B
Uniprot ID	CDY2_HUMAN
Immunogen	
Immunogen Region	1-280 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CDY2B (NP_001001722.1).
Immunogen Sequence	MASQEFVEAIVDKRQDKNG NTQYLVRWKGYDKQDDTWEP EQHLMNCEKCVHDFNRRQTE KQKCLTWTTTSRIFSNNARR RTSRSTKANYSKNSPKTPVT DKHHRSKNCKLFAASKNVRR KAASTLSDTKNMEIINSTIE TLAPDSPFDHKKTVSGFQKL EKLDPAAADQQDTVVKVTE GKLLRDPLSHPGAEGTGIQN KTQMHPLMSQMSGSVTASMA TGSATRKGVVLIDPLAAN