

Anti-ZC3H12A antibody (1-140 aa) (STJ11106072)

STJ11106072

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	80149
Gene Symbol	ZC3H12A
Uniprot ID	ZC12A_HUMAN
Immunogen	
Immunogen Region	1-140 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZC3H12A (NP_079355.2).
Immunogen Sequence	MSGPCGKPKVLEASPTMSLW EFEDSHSRQGTTPRPGQELAA EEASALELQMKVDFFRKLGY SSTEIHSVLQKLGVQADTNT VLGLVVKHGTATERERQTSP DPCPQLPLVPRGGGTPKAPN LEPPLPEEEKEGSDLRPVVI