

Anti-SCG3 antibody (199-468 aa) (STJ11106075)

STJ11106075

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	29106
Gene Symbol	SCG3
Uniprot ID	SCG3_HUMAN
Immunogen	
Immunogen Region	199-468 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 199-468 of human SCG3 (NP_037375.2).
Immunogen Sequence	KLISKEANNYEEDPNKPTSW TENQAGKIPEKVTMAAIQD GLAKGENDETVSNTLTNG LERRTKTYSEDNFEELQYFP NFYALLKSIDSEKEAKEKET LITIMKTLIDFVKMMVKYGT ISPEEGVSYLENLDEMIALQ TKNKLEKNATDNISKLPAP SEKSHEETDSTKEEAAKMEK EYGLKDDSTKDDNSNPGGKT DEPKGKTEAYLEAIRKNIEW LKKHDKKGNKEDYDLKMR