

Anti-IFN-beta antibody (22-187 aa) (STJ11106080)

STJ11106080

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3456
Gene Symbol	IFNB1
Uniprot ID	IFNB_HUMAN
Immunogen	
Immunogen Region	22-187 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-187 of human IFN-beta (NP_002167.1).
Immunogen Sequence	LQRSSNFQCQKLLWQLNGRL EYCLKDRMNFDIPEEIKQLQ QFQKEDAALTIYEMLQNIFA IFRQDSSSTGWNENETIVENLL ANVYHQINHLKTVLEEKLEK EDFTRGKLMSSLHLKRYYYGR ILHYLKAKEYSHCAWTIVRV EILRNFYFINRLTGYLRLN