

Anti-KIAA1456 antibody (110-330 aa) (STJ11106086)

STJ11106086

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	57604
Gene Symbol	TRMT9B
Uniprot ID	TRM9B_HUMAN
Immunogen	
Immunogen Region	110-330 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 110-330 of human KIAA1456 (NP_065895.2).
Immunogen Sequence	VIHHFSTKQRRIRAIKEMAR VLVPGGQLMIYVWAMEQKNR HFEKQDVLVPWNRALCSQLF SESSQSGRKRCQCGYPERGHP YHPPCSECSVCVFKEQCGS KRSHSVGYEPAMARTCFANI SKEGEEYGFYSTLGKSFRS WFFSRSLDESTLRKQIERVR PLKNTVEVWASSTVTVPQSRH SSLDFDHQEPFSTKGQSLDE EVFVSSSGKHLEWLAPGT L