

Anti-CPB2 antibody (115-365 aa) (STJ11106088)

STJ11106088

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1361
Gene Symbol	CPB2
Uniprot ID	CBPB2_HUMAN
Immunogen	
Immunogen Region	115-365 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 115-365 of human CPB2 (NP_001863.2).
Immunogen Sequence	ASASYEQYHSLNEIYSWIE FITERHPDMLTKIHIGSSF KYPLYVLKVGKEQAAKNAI WIDCGIHAREWISPAFCLWF IGHITQFYGIIGQYTNLLRL VDFYVMPVNVNVDGYDSWKK NRMWRKNRSFYANNHCIGTD LNRNFASKHWCEEGASSSC SETYCGLYPESEPEVKAVAS FLRRNINQIKAYISMHSYSQ HIVFPYSYTRSKSKDHEELS LVASEAVRAIEKTSKNTRY