

Anti-SEMA5B antibody (800-980 aa) (STJ11106093)

STJ11106093

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	54437
Gene Symbol	SEMA5B
Uniprot ID	SEM5B_HUMAN
Immunogen	
Immunogen Region	800-980 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 800-980 of human SEMA5B (NP_001243276.1).
Immunogen Sequence	RRACENGNSCLGCGVEFKTC NPEGCPVRRNTPWTPWLPV NVTQGGARQEQRRFTCRAP LADPHGLQFGRRTETRTCP ADGSGSCDTDALVEVLLRSG STSPHTVSGGWAAGPWSSC SRDCELGFRVRKRTCTNPEP RGGGLPCVGDAAEYQDCNPQ ACPVRGAWSCWTSWSPCSAS C