

Anti-F8 antibody (335-630 aa) (STJ11106099)

STJ11106099

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2157
Gene Symbol	F8
Uniprot ID	FA8_HUMAN
Immunogen	
Immunogen Region	335-630 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 335-630 of human F8 (NP_000123.1).
Immunogen Sequence	QHDGMEAYVKVDSCPEEPQL RMKNNEEAEDYDDDLTDSEM DVVRFDDDNSPSFIQIRVA KKHPKTWVHYIAEEEEEDWDY APLVLAPDDRSYKSQYLNNG PQRIGRKYKKVRFMAYTDET FKTREAIQHESGILGPLLYG EVGDTLLIIFKNQASRPYNI YPHGITDVRPLYSRRLPGKV KHLKDFPILPGEIFKYKWTV TVEDGPTKSDPRCLTRYSS FVNMERDLASGLIGPLLIC