

Anti-TCOF1 antibody (659-958 aa) (STJ11106106)

STJ11106106

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6949
Gene Symbol	TCOF1
Uniprot ID	TCOF_HUMAN
Immunogen	
Immunogen Region	659-958 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 659-958 of human TCOF1 (NP_001008657.1).
Immunogen Sequence	VRVGTQAPRKAGTATSPAGS SPAVAGGTQRPAAEDSSSSEE SDSEEEKTGLAVTVGQAKSV GKGLQVKAASVPVKGSLGQG TAPVLPGKTGPTVTQVKA EK QEDSESSEESDSEEAASP AQVKTSVKKTQAKANPAAAR APSAKGTISAPGKVVTAQAQ AKQRSPSKVKPPVRNPQNST VLARGPASVPSVGKAVATAA QAQTGPEEDSGSSEESDSE EEAETLAQVKPSGKTHQIR