

Anti-p60 CAF1 antibody (270-559 aa) (STJ11106107)
STJ11106107

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	8208
Gene Symbol	CHAF1B
Uniprot ID	CAF1B_HUMAN
Immunogen	
Immunogen Region	270-559 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 270-559 of human p60 CAF1 (NP_005432.1).
Immunogen Sequence	SRKNLKRPIAHLPCPGKATL AVRCCPVYFELRPVETGVE LMSLPYRLVFAVASEDSVLL YDTQQSFPGYVSNIIHYHTL SDISWSSDGAFLAISSTDGY CSFVTFEKDELGIPLKEKPV LNMRTPDATKTKSQTHRGS SPGPRPVEGTPASRTQDPSS PGTTPQARQAPAPTIVIRDP PSITPAVKSPPLGPGSEKTL QPSSQNTKAHPSRRVTLNLT QAWSKTTPRRINLTPLKTD