

Anti-MED23 antibody (1046-1376 aa) (STJ11106109)

STJ11106109

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	70208
Gene Symbol	Med23
Uniprot ID	MED23_MOUSE
Immunogen	
Immunogen Region	1046-1376 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1046-1376 of mouse MED23 (NP_001159888.1).
Immunogen Sequence	NRPGWCLSDTYLKHAMNAR EDNPWVPEDSYCKLIGRLV DTMAGKSPGPFPCDWRFE FPNPAAHALHVTCELMALA VPGKDVGNALLNVVLKSQPL VPRENITAWMNAIGLIITALEPEYWIVLHDRVNVISSSS LTSETEWVGYPFRFLDFTAC HQSYSSEMSCSYTLALAHAVW HHSSIGQLSLIPKFLTEALL PVVKTEFQLLYVYHLVGPFL QRFQQRTRCMIEIGVAFY