

Anti-PLCD4 antibody (1-290 aa) (STJ11106116)

STJ11106116

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84812
Gene Symbol	PLCD4
Uniprot ID	PLCD4_HUMAN
Immunogen	
Immunogen Region	1-290 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PLCD4 (NP_116115.1).
Immunogen Sequence	MASLLQDQLTTDQDLLLMQE GMPMRKVRSKSWKKLRYFRL QNDGMTVWHARQARGSAKPS FSISDVETIRNGHDSSELLRS LAELPLEQGFTIVFHGRRS NLDLMANSVEEAQIWMRGLQ LLVDLVTSMDHQERLDQWLS DWFQRGDKNQDGKMSFQEVQ RLHLMLNVEMDQEYAFSLFQ AADTSQSGTLEGEFVQFYK ALTKRAEVQELFESFSADGQ KLTLLFLDFLQEEQKERD