

Anti-CXCL3 antibody (1-107 aa) (STJ11106120)

STJ11106120

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2921
Gene Symbol	CXCL3
Uniprot ID	CXCL3_HUMAN
Immunogen	
Immunogen Region	1-107 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human CXCL3 (NP_002081.2).
Immunogen Sequence	MAHATLSAAPSNPRLLRVAL LLLLLVAASRRRAAGASVVTE LRCQCLQTLQGIHLKNIQSV NVRSPGPHCAQTEVIATLKN GKKACLNPASPMVQKIIIEKI LNKGSTN