

Anti-NLGN2 antibody (15-180 aa) (STJ11106122)

STJ11106122

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	57555
Gene Symbol	NLGN2
Uniprot ID	NLGN2_HUMAN
Immunogen	
Immunogen Region	15-180 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 15-180 of human NLGN2 (NP_065846.1).
Immunogen Sequence	QRGGGGPGGGAPGGPGLGLG SLGEERFPVVNTAYGRVRGV RRELNNELGPVQFLGVPI ATPPLGARRFQPPEAPASWP GVRNATTLPPACPQNLHGAL PAIMLPWWFTDNLEAAATYV QNQSEDCLYLNLYVPTEDGP LTKKRDEATLNPPDTPDIRDP GKPKVM