

Anti-RPUSD3 antibody (32-351 aa) (STJ11106124)

STJ11106124

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	285367
Gene Symbol	RPUSD3
Uniprot ID	RUSD3_HUMAN
Immunogen	
Immunogen Region	32-351 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 32-351 of human RPUSD3 (NP_775930.2).
Immunogen Sequence	PEDAGFGTEARHQRQPRGSC QRSGLDQPFAGLLPKNLS REELVDALRAAVVDRKGPLV TLNKPQGLPVTGKPGELTLF SVLPESQSLGLREQELQVV RASGKESSGLVLLSSCPQTA SRLQKYFTHARRAQRPTATY CAVTDGIPAASEGKIQAALK LEHIDGVNLTVPVKAPSRKD ILEGVKKTLSHFRVVGSG CALVQLPLTVFSSQLQVHM VLQLCPVLGDHMY SARVGT