

Anti-NOV antibody (71-357 aa) (STJ11106126)

STJ11106126

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4856
Gene Symbol	CCN3
Uniprot ID	CCN3_HUMAN
Immunogen	
Immunogen Region	71-357 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 71-357 of human NOV (NP_002505.1).
Immunogen Sequence	RGESCDLEPCDESSGLYCD RSADPSNQTGICTAVEGDNC VFDGVIYRSGEKFQPSCKFQ CTCRDGQIGCVPRCQLDVLL PEPNCAPARKVEVPGECCCK WICGPDEEDSLGGLTLAAYR PEATLGVEVSDSSVNCIEQT TEWTACSKSCGMGFSTRVTN RNRQCEMLKQTRLCMVRPCE QEPEQPTDKKGKKCLRTKKS LKAIHLQFKNCTSLHTYKPR FCGVCSDGRCCTPHNTKTI