

Anti-TLR3 antibody (740-904 aa) (STJ11106127)

STJ11106127

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	
Immunogen Region	740-904 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 740-904 of human TLR3 (NP_003256.1).
Immunogen Sequence	RVLGFKEIDRQTEQFEYAAY IIHAYKDKDWVWEHFSSMEK EDQSLKFCLEERDFEAGVFE LEAIVNSIKRSRKIIIFVITH HLLKDPLCKRFKVVHHAQQAA IEQNLDIILVFLEIIPDYK LNHALCLRRGMFKSHCILNW PVQKERIGAFRHKLQVALGS KNSVH