

Anti-ZDHHC14 antibody (273-488 aa) (STJ11106130)
STJ11106130

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	79683
Gene Symbol	ZDHHC14
Uniprot ID	ZDH14_HUMAN
Immunogen	
Immunogen Region	273-488 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 273-488 of human ZDHHC14 (NP_078906.2).
Immunogen Sequence	TYLISSNQTTNEDIKGSWSN KRGKENYNPYSGNIFTNCC VALCGPISPSLIDRRGYIQP DTPQPAAPSNGITMYGATQS QSDMCDQDQCISTKFVLQA AATPLLQSEPSLTDELHLP GKPGLGTPCASLTLPPTPP ASMPNLAEATLADVMPRKDE HMGHQFLTPEAPSPRLLA AGSPLAHSRTMHVLGLASQD SLHEDSVRGLVKLSSV